1

2

3

4

FR

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UK

Qflow® sensors are single-use sensors specifically developed to operate with the SPIROLYSER® Q13® spirometer. Qflow® sensors are designed and adapted for lung function tests (spirometry). Qflow® sensors are designed so that their combined resistance to the spirometer Spirolyser® Q13® does not exceed the ATS / ERS recommendations.

The technologically innovative design of the Spirolyser® Q13® and the Qflow® sensor minimise the risk of cross-contamination between patients. Bacteriological tests performed by the **Public Health England** (Salisbury - UK) laboratory show that the internal design of the product, associated with the Qflow® sensor, decreases the risk of patient cross-contamination to 99.999% if the Qflow® sensor is changed between each patient. Its innovative design also avoids contamination in the non-accessible internal parts.

Position the Single-Use Qflow® Sensor:

Step 1, 2, 3

- Locate the dotted line opening part
- Hold everything with one hand by the mouthpiece side (small side of packaging)
- With the other hand, open the detachable part of the Qflow® packaging along the dotted line
- Insert the Qflow® into the opening of the sensor until it stops
- Remove and discard the rest of the packaging

Eject the Single-Use Qflow® Sensor:

Step 4

- Turn the device downwards, above a rubbish bin
- Push the trigger upwards, located on the back of the spirometer handle
- The Qflow® sensor is ejected

Caution: Qflow® sensors are single use.

Caution: The use of single-use sensors does not exempt from decontaminating the equipment after each use.

Caution: Used or dirty accessories must be disposed of with biological waste (DASRI).

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