
Garantie  
sans Latex

Usage unique

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**UK HEALTH SECURITY AGENCY**  
Porton Down  
Salisbury  
United Kingdom

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Guaranteed  
Latex free

Single use

The Ergofilter® SP1 (with mouthpiece) and SP1M (with flexible jaw tip) are single-use antibacterial and antiviral filters, intended to be used with

**Caution:**

- Always use an *Ergofilter® SP1* or *SP1M* filter between the patient and the spirometer.
- *Ergofilter® SP1* and *SP1M* filters are single-use.
- The use of filters does not exempt Fim médical spirometer from regular disinfection.
- Used or dirty accessories must be disposed of with biological waste (DASRI) to avoid any risk of cross-contamination.
- Any serious incident in connection with the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

the spirometers of the FIM MEDICAL brand.

Ergofilter® filters are designed and adapted for respiratory functional exploration exams (spirometry). Ergofilter® filters are designed so that their cumulative resistance to that of the spirometer does not exceed the recommendations of the ATS / ERS.

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CE marking Directive 93/42/EEC or MDR 2017/745 + N° of notified body

Consult the accompanying documents

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<sup>2</sup> The bacterial and viral efficacy tests of Ergofilter® filters were carried out by the laboratory:

UK HEALTH SECURITY AGENCY  
Porton Down  
Salisbury  
United Kingdom

**Installation of the Ergofilter® SP1 or SP1M filters:**

- Take the filter out of its unit packaging
- Do not touch the mouth of the filter with your fingers but hold it by the circumference
- Use the filter in the right direction because it is not reversible
- Mount the cone on the spirometer (label side "tip") then connect the filter to the cone

**Planned operators:**

The SP1 and SP1M antibacterial filters should only be fitted by operators of the spirometer. The operators are health professionals such as medical secretaries, nurses, doctors, pulmonologists, or any other specialist doctors who have received recognized and comprehensive training in the use of spirometers.

**Patient population**

The patient population that can use this device is:

- Persons aged 4 - 95 years.
- Wishing to perform spirometry with an Fim médical spirometer to assess their breathing capacities.
- All people not suffering from a contraindication related to their respiratory condition

**Indication: Antibacterial and antiviral filters to protect against bacterial and viral risks between patients during spirometry.**

**Medical contraindications to performing spirometry:**

- Current or recent pneumothorax
- Recent pleural puncture or biopsy
- Hemoptysis in progress
- Severe asthma attack
- Bacilliferous tuberculosis (risk to staff)
- Bronchial infections, infectious lung disease
- Acute asthma attack
- Decompensation of chronic respiratory failure
- Inability to perform breathing maneuvers; ventilatory motor incoordination
- Abnormal pain and fatigue
- Myocardial infarction <6 months
- Laparotomy <6 months